

EVI 91 Water Bath



Inlabo water bath for rotary evaporator EVI

GENERAL CHARACTERISTICS

- ❑ Built in stainless steel AISI 304 with a double side both for a better thermic isolation from the heating element and to protect from scorches.
- ❑ Bath completely free inside and easy to clean.
- ❑ Bath diameter 240 mm, stainless steel AISI 304/L (18/10).
- ❑ Capacity 4,3 litres.
- ❑ Suitable for flasks up to 3 litres.
- ❑ Power requirements: 220V 50 HZ 1500 Watt.
- ❑ ON/OFF thermostat for electronic temperature regulation.
- ❑ Regulation accuracy: $\pm 1,5$ °C.
- ❑ Control range 20 °C – 95 °C.
- ❑ Display: - 50 °C + 150 °C.
- ❑ Safety: limit temperature 100 °C (+0/-6 °K) with limit thermostat, manual reset.
- ❑ Size: mm 288x276x288.
- ❑ Weight: Kg 8.

